

Code No.05-0437

SPECIFICATION

o-Ethoxybenzaldehyde

Wako Special Grade

| REQUIREMENT | SPECIFICATION |
|---|---|
| Appearance | Colorless ~ slightly yellow, clear liquid |
| Density (20°C) | 1.081~1.086g/ml |
| Refractive index 20°C | 1.541~1.545 |
| Water | max.0.1% |
| Acidity (as C ₆ H ₅ COOH) | max.0.5% |
| Assay (capillary GC) | min.97.0% |