

SPECIFICATION

Diethyl Ether for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 300

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay (C ₂ H ₅ OC ₂ H ₅) (GC)	min.99.5%(mass/mass)
Density (20°C)	0.712~0.714g/ml
Water	max.0.1%(mass/mass)
Residue after evaporation	max.0.0002%(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test
* Suitability for pesticide residue analysis (GC-FPD)	to pass test

JIS K 8357:2012

*:Additional test performed by Wako but not required by JIS