

## SPECIFICATION

## Ethyl Acetate

for Preparative Chromatography

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.897~0.902g/ml
Refractive index 20°C	1.372~1.374
Absorbance (260nm)	max.0.40
Absorbance (300~400nm)	max.0.01
Water	max.0.1%
Residue after evaporation	max.5ppm
Acidity (as CH <sub>3</sub> COOH)	max.0.003%
Assay (capillary GC)	min.99.7%