

Code No.05-0560

## SPECIFICATION

## 2-Ethylhexyl Methacrylate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.881~0.887g/mL
Refractive index 20°C	1.437~1.443
Water	max.0.1%
Acidity [as CH <sub>2</sub> :C(CH <sub>3</sub> )COOH]	max.0.1%
Assay (capillary GC)	min.98.0%