

Code No.05-0618

SPECIFICATION

Ethyl Acetate, Dehydrated
for Organic Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.898~0.902g/ml
Refractive index 20°C	1.370~1.374
Water	max.0.005%
Residue after evaporation	max.0.002%
Acidity (as CH ₃ COOH)	max.0.005%
Aldehydes	試験適合
Substances darkened by sulfuric acid	to pass test
Assay (capillary GC)	min.99.5%