

SPECIFICATION

N-Ethyldiisopropylamine

for Peptide Synthesis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Water	max.0.01%
Ethylamine (CH ₃ CH ₂ NH ₂)	max.0.0015%
N-Ethylisopropylamine [(CH ₃) ₂ CHNHC ₂ H ₅]	max.0.015%
Isopropylamine [(CH ₃) ₂ CHNH ₂]	max.0.005%
Assay (capillary GC)	min.99.0%