

Code No.05-0637

SPECIFICATION

86% Ethanol-ME, Denaturated

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.785~0.791g/ml
Refractive index 20°C	1.353~1.359
Water	max.0.2%
Acidity (as CH ₃ COOH)	max.0.002%
Substances darkened by sulfuric acid	to pass test
Ethanol	85.0~87.0%
Methanol	12.5~14.0%
2-Propanol	0.5~1.5%