

## SPECIFICATION

## Ethyl Acetate for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
** Assay (CH <sub>3</sub> COOC <sub>2</sub> H <sub>5</sub> ) (GC)	min.99.7%(mass/mass)
Density (20°C)	0.899~0.902g/mL
Water	max.0.05%(mass/mass)
Residue after evaporation	max.0.0002%(mass/mass)
Suitability for pesticide residue-PCB analysis	to pass test
* Suitability for pesticide residue analysis (GC-FPD)	to pass test

JIS K 8110:2012

\*:Additional test performed by Wako but not required by JIS

\*\*:The specification is originally set by Wako not JIS