

SPECIFICATION

Ethyl Acetate for Pesticide Residue And Polychlorinated Biphenyl Analysis
Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance ** Assay (CH ₃ COOC ₂ H ₅) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis * Suitability for pesticide residue analysis (GC-FPD)	Colorless clear liquid min.99.7%(mass/mass) 0.899~0.902g/mL max.0.05%(mass/mass) max.0.0002%(mass/mass) to pass test to pass test

JIS K 8110:2012

*:Additional test performed by Wako but not required by JIS

**:The specification is originally set by Wako not JIS