

SPECIFICATION

Ethanol(99.5) for Pesticide Residue And Polychlorinated Biphenyl Analysis

Concentration 5000

REQUIREMENT	SPECIFICATION
* Appearance Assay (C ₂ H ₅ OH) (GC) Density (20°C) Water Residue after evaporation Suitability for pesticide residue-PCB analysis * Suitability for pesticide residue analysis (GC-FPD)	Colorless clear liquid min.99.5%(mass/mass) 0.789~0.791g/mL max.0.3%(mass/mass) max.0.0002%(mass/mass) to pass test to pass test

JIS K 8093:2012

*:Additional test performed by Wako but not required by JIS