

SPECIFICATION

Ethanol (99.5)
for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Density (20°C)	0.789~0.791g/mL
Water	max.0.4%
Acidity (as CH ₃ COOH)	max.0.002%
Alkalinity (as NH ₃)	max.1 ppm
Substances darkened by sulfuric acid	to pass test
Deoxyribonuclease	to pass test
Ribonuclease	to pass test
Assay (capillary GC)	min.99.5%