

SPECIFICATION

Formic Acid

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Assay (HCOOH) Solubility in water Residue after evaporation Acetic acid (CH ₃ COOH) (GC) Chloride (Cl) Sulfate (SO ₄) Sulfite Copper (Cu) Lead (Pb) Iron (Fe)	Colorless clear liquid min.98.0%(mass/mass) to pass test max.0.002%(mass/mass) max.0.3%(mass/mass) max.5ppm(mass/mass) max.0.005%(mass/mass) to pass test max.3ppm(mass/mass) max.3ppm(mass/mass) max.3ppm(mass/mass)

JIS K 8264:2020(*:Additional test performed by Wako but not required by JIS)The assay of this product is approximately 98.0%.