

# SPECIFICATION

## Formic Acid (abt.99%) for LC/MS

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Solubility in water	to pass test
Absorbance (1→4, 254nm)	max.1.00
Absorbance (1→4, 260nm)	max.0.20
Absorbance (1→4, 270nm)	max.0.02
Absorbance (1→4, 300~400nm)	max.0.01
Water	max.1.0%
Fluorescence test	to pass test
Suitability for LC/MS analysis	to pass test
Assay (HPLC)	min.99.5%