

Code No.06-0525

SPECIFICATION

Florisol, 150~250 μ m (60~100mesh)
for Column Chromatography

| REQUIREMENT | SPECIFICATION |
|--|---------------------|
| Appearance | to pass test |
| Particle size (150~250 μ m) | min.60% |
| Particle size (coarser than 250 μ m) | max.40% |
| Particle size (finer than 150 μ m) | max.20% |
| Bulk density | under investigation |
| Loss on drying at 180°C | max.2.0% |
| pH (50g/l aqueous slurry,25°C) | 8.6~10.2 |
| Methanol dipping test | to pass test |