

SPECIFICATION

Florisol, 75~150 μ m (100~200mesh)

for Column Chromatography

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Particle size (75~150 μ m)	min.60%
Particle size (coarser than 150 μ m)	max.20%
Particle size (finer than 75 μ m)	max.40%
Bulk density	0.4~0.7g/ml
Loss on drying at 180°C	max.2.0%
pH (50g/l aqueous slurry,25°C)	8.6~10.2
Methanol dipping test	to pass test