

Code No.06-0633

# SPECIFICATION

## Fosfomycin Sodium Standard for High Performance Liquid Chromatography

REQUIREMENT	SPECIFICATION
Appearance	White ~ slightly brown, crystals ~ powder
Solubility in water	to pass test
Specific rotation [ $\alpha$ ]D <sub>20</sub> (c=5,H <sub>2</sub> O) (calculated on the dehydrous basis)	-3.5~-5.5°
Water	max.3.0%
Assay (HPLC)	min.95.0%