

Code No.06-0653

## SPECIFICATION

Fumonisin Mixture Standard Solution [B1, B2 each 50  $\mu\text{g/mL}$  Acetonitrile:Water(1:1)  
Solution]  
for Mycotoxin Assay

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Practical test	to pass test
Fumonisin B1	49.2~50.8 $\mu\text{g/mL}$
Fumonisin B2	49.2~50.8 $\mu\text{g/mL}$