

Code No.07-0501

SPECIFICATION

D(+)-Galactosamine Hydrochloride
for Molecular Biology

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, powder
Solubility in isotonic sodium chloride solution	to pass test
Specific rotation [α]D20(c=1,H2O,24hr)	+90~+100°
Water	max.1.0%
Residue after ignition (as sulfate)	max.1.0%
Deoxyribonuclease	to pass test
Ribonuclease	to pass test
Thin layer chromatography test	to pass test