

Code No.07-0523

SPECIFICATION

Ginsenoside Rc

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description (JP method)	White crystalline powder
Purity (JP method)	to pass test
* Solubility in methanol	to pass test
* Loss on drying over silica gel in vacuum	max.2.0%
* Assay (HPLC)	min.98.0%

This product is suited for Ginsenoside Rc of general tests on The Japanese Pharmacopoeia.(*:Additional test performed by Wako)