

Code No.07-0633

SPECIFICATION

(+)-Galocatechin Standard  
for Food Analysis

REQUIREMENT	SPECIFICATION
Appearance	White ~ pale brown,
Solubility in methanol	crystalline powder ~ powder
	to pass test
Specific rotation [ $\alpha$ ]D20(c=0.2,CH3OH)	under investigation
(calculated on the dehydrous basis)	
Water	max.3.0%
Assay (HPLC)	min.99.0%
Assay (qNMR) (as galocatechin)	min.90.0%