

Code No.07-0670

## SPECIFICATION

 $\alpha$ -Glucosylrutin

for Biochemistry

REQUIREMENT	SPECIFICATION
Appearance	Pale yellow ~ yellow, powder
Solubility in water	to pass test
pH of a 10g/L solution at 25°C	4.5~6.5
Loss on drying at 120°C	max.6.0%
Residue after ignition (as sulfate)	max.2.0%
Quercetin	under investigation
Assay (as rutin) (calculated on the dried basis)	min.70.0%