

Code No.08-0015

## SPECIFICATION

## Heptane

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Colorless clear liquid
Assay [CH <sub>3</sub> (CH <sub>2</sub> ) <sub>5</sub> CH <sub>3</sub> ] (GC)	min.99.0%(mass/mass)
Density (20°C)	0.682~0.686g/ml
Refractive index 20°C	1.386~1.389
Residue after evaporation	max.0.001%(mass/mass)
Substances darkened by sulfuric acid	to pass test

JIS K 9701:2013(\*:Additional test performed by Wako but not required by JIS)