

Code No.08-0105

SPECIFICATION

Hydrobromic Acid

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly brown, clear liquid
Chloride (Cl)	max.0.25%
Phosphate (PO ₄)	max.0.002%
Sulfate (SO ₄)	max.0.02%
Sulfite	to pass test(SO ₃ :max.about0.006%)
Heavy metals (as Pb)	max.0.002%
Iron (Fe)	max.0.001%
Assay	47.0~49.0%