

# SPECIFICATION

## Hydrochloric Acid for Analysis of Poisonous Metals

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Zinc (Zn)	max.2ppb
Cadmium (Cd)	max.0.5ppb
Mercury (Hg)	max.0.1ppb
Tin (Sn)	max.2ppb
Lead (Pb)	max.0.5ppb
Arsenic (As)	max.2ppb
Antimony (Sb)	max.1ppb
Vanadium (V)	max.10ppb
Selenium (Se)	max.2ppb
Chromium (Cr)	max.0.5ppb
Manganese (Mn)	max.2ppb
Assay	35.0~37.0%