

# SPECIFICATION

## Hydrogen Peroxide

for Atomic Absorption Spectrochemical Analysis

REQUIREMENT	SPECIFICATION
Appearance	Colorless clear liquid
Copper (Cu)	max.1ppb
Silver (Ag)	max.1ppb
Magnesium (Mg)	max.10ppb
Calcium (Ca)	max.50ppb
Zinc (Zn)	max.10ppb
Cadmium (Cd)	max.1ppb
Mercury (Hg)	max.1ppb
Tin (Sn)	max.10ppb
Lead (Pb)	max.5ppb
Arsenic (As)	max.1ppb
Chromium (Cr)	max.1ppb
Manganese (Mn)	max.5ppb
Iron (Fe)	max.5ppb
Assay	30.0~35.5%