

Code No.08-0438

## SPECIFICATION

## 2-Hydroxyethyl Methacrylate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	1.069~1.077g/ml
Refractive index 20°C	1.450~1.456
Water	max.0.5%
Acidity [as CH <sub>2</sub> :C(CH <sub>3</sub> )COOH]	max.1.0%
Assay (capillary GC)	min.95.0%