

Code No.08-0481

## SPECIFICATION

## Heptaethylene Glycol Monododecyl Ether

Wako Special Grade

| REQUIREMENT                        | SPECIFICATION   |
|------------------------------------|-----------------|
| Appearance                         | to pass test    |
| Solubility in ethanol              | to pass test    |
| Density (35°C)                     | 0.977~0.983g/ml |
| Refractive index n <sub>D</sub> 35 | 1.448~1.454     |
| Water                              | max.0.5%        |
| Assay (capillary GC)               | min.98.0%       |