

Code No.08-0607

SPECIFICATION

Heptyl Acetate

Wako 1st Grade

| REQUIREMENT | SPECIFICATION |
|-----------------------------------|--|
| Appearance | Colorless ~ slightly yellow, clear liquid |
| Density (20°C) | 0.866~0.874g/mL |
| Refractive index 20°C | 1.412~1.417 |
| Water | max.0.2% |
| Acidity (as CH ₃ COOH) | max.0.1% |
| Assay (capillary GC) | min.98.0% |