

Code No.08-0846

SPECIFICATION

Hirsutine

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description (JP method)	White ~ pale yellow, crystals or crystalline powder
Optical rotation [α]D20 (JP method)	about +70°
Absorbance E1%1cm (245nm) (JP method)	354~379
Purity (JP method)	-
Analogous substances (JP method)	to pass test
Solubility in methanol	to pass test
Water	under investigation
Assay (HPLC)	min.98.0%