

Code No.08-0869

# SPECIFICATION

## (E)-Isoferulic Acid

for the Japanese Pharmacopoeia General Tests

REQUIREMENT	SPECIFICATION
Description (JP method)	White to light yellowcrystals or crystalline powder
Identification (JP method)	to pass test
Purity (JP method)	-
Related substances (JP method)	to pass test
* Solubility in methanol	to pass test
* Water	measured value
* Assay (HPLC)	min.98.0%

This product is suited for (E)-Isoferulic Acid on The Japanese Pharmacopoeia.(\*:Additional test performed by Wako)