

Code No.08-1045

# SPECIFICATION

## Hyaluronic Acid Sodium Salt, from Rooster Comb for Biochemistry

REQUIREMENT	SPECIFICATION	
Appearance	White ~ pale yellow, crystalline powder ~ powder or granules to pass test	
Solubility in water		
Loss on drying at 60°C over phosphorus(V) oxide in vacuum		max.15.0%
Residue after ignition (as sulfate) (after drying)		15.0~20.0%
pH of a 1g/L solution at 25°C		5.0~8.5
Nitrogen compounds (as N)		2.8~4.0%