

Code No.08-1076

SPECIFICATION

L-Histidine

Japanese Pharmacopoeia

REQUIREMENT	SPECIFICATION
* Description (JP)	White crystals or crystalline powder
Identification (JP)	to pass test
Optical rotation [α]D20 (JP)	+11.8~+12.8°
pH (JP)	7.0~8.5
Purity (JP)	-
(1) Clarity and color of solution (JP)	to pass test
(2) Chloride (JP)	not more than 0.021%
(3) Sulfate (JP)	not more than 0.028%
(4) Ammonium (JP)	not more than 0.02%
(5) Iron (JP)	not more than 10ppm
(6) Related substances (JP)	to pass test
Loss on drying (JP)	not more than 0.3%
Residue on ignition (JP)	not more than 0.1%
Assay (calculated on the dried basis) (JP)	99.0~101.0%
* Bacterial endotoxins	less than 2.0EU/g
Identification (USP-NF)	to pass test
Assay (dried basis) (USP-NF)	98.5~101.5%
Impurities (USP-NF)	-
(1) Residue on ignition (USP-NF)	NMT 0.4%
(2) Chloride (USP-NF)	NMT 0.05%
(3) Sulfate (USP-NF)	NMT 0.03%
(4) Iron (USP-NF)	NMT 30ppm
(5) Related compounds (USP-NF)	to pass test
Specific tests (USP-NF)	-
(1) Optical rotation (USP-NF)	+12.6~+14.0°
(2) pH (USP-NF)	7.0~8.5
(3) Loss on drying (USP-NF)	NMT 0.2%
Identification (Ph.Eur.)	to pass test
Appearance of solution (Ph.Eur.)	to pass test

Specific optical rotation (Ph.Eur.)	+11.4~+12.4°
Ninhydrin-positive substances (Ph.Eur.)	-
(1) For each impurity (Ph.Eur.)	maximum 0.2%
(2) Total (Ph.Eur.)	maximum 0.5%
Chlorides (Ph.Eur.)	maximum 200ppm
Sulfates (Ph.Eur.)	maximum 300ppm
Ammonium (Ph.Eur.)	maximum 0.02%
Iron (Ph.Eur.)	maximum 10ppm
Loss on drying (Ph.Eur.)	maximum 0.5%
Sulfated ash (Ph.Eur.)	maximum 0.1%
Assay (dried substance) (Ph.Eur.)	98.5~101.0%

Japanese Pharmacopoeia (L-Histidine)(*:Additional test performed by Wako)United States Pharmacopeia - National Formulary (Histidine)European Pharmacopoeia (Histidine)