

Code No.09-0457

## SPECIFICATION

## Indium(III) Nitrate Trihydrate

REQUIREMENT	SPECIFICATION
Appearance	White ~ nearly white, crystals ~ powder or mass to pass test
Solubility in dil.hydrochloric acid	
pH of a 50g/L solution at 25°C	2.0~3.0
Chloride (Cl)	max.0.001%
Copper (Cu)	max.0.001%
Lead (Pb)	max.0.001%
Iron (Fe)	max.0.001%
Assay	min.97.0%