

Code No.09-0470

# SPECIFICATION

## Interferon- $\gamma$ , Mouse, recombinant for Biochemistry

| REQUIREMENT                   | SPECIFICATION             |
|-------------------------------|---------------------------|
| Appearance                    | to pass test              |
| Bacterial endotoxins          | less than 0.1 ng/ $\mu$ g |
| Electrophoresis test          | to pass test              |
| Mass spectrometry (MALDI-TOF) | to pass test              |

Source: E. coli expressed mouse interferon- $\gamma$