

Code No.09-0515

## SPECIFICATION

## Isopropyl Myristate

Wako 1st Grade

REQUIREMENT	SPECIFICATION
Appearance	Colorless ~ slightly yellow, clear liquid
Density (20°C)	0.850~0.857g/ml
Refractive index 20°C	1.430~1.437
Water	max.0.1%
Acidity [as CH <sub>3</sub> (CH <sub>2</sub> ) <sub>12</sub> COOH]	max.0.2%
Assay (capillary GC)	min.95.0%