

Code No.09-0613

SPECIFICATION

Interleukin-3, Mouse, recombinant, Animal-derived-free  
for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01 ng/ μg
Electrophoresis test (SDS-PAGE)	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source:E. coli expressed mouse interleukin-3  
This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.