

Code No.09-0698

SPECIFICATION

Interferon- γ , Mouse, recombinant, Animal-derived-free
for Cellbiology

REQUIREMENT	SPECIFICATION
Appearance	to pass test
Bacterial endotoxins	less than 0.01ng/ μ g
Electrophoresis test (SDS-PAGE)	to pass test
Mass spectrometry (MALDI-TOF)	to pass test

Source:E. coli expressed mouse interferon- γ
This product is recombinant cytokine which is produced in the process of culture and purification without using any substances from animals.