

Code No.12-0037

SPECIFICATION

Lead

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance Copper (Cu) Silver (Ag) Zinc (Zn) Tin (Sn) Arsenic (As) Bismuth (Bi) Iron (Fe)	Silver gray ~ bluish gray, shot max.5ppm max.0.001% max.0.001% max.0.001% max.1ppm max.0.002% max.0.001%

JIS K 8701-1994

*:Additional test performed by Wako but not required by JIS
This product is shot