

Code No.12-0037

SPECIFICATION

Lead

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	Silver gray ~ bluish gray, shot
Copper (Cu)	max.5ppm
Silver (Ag)	max.0.001%
Zinc (Zn)	max.0.001%
Tin (Sn)	max.0.001%
Arsenic (As)	max.1ppm
Bismuth (Bi)	max.0.002%
Iron (Fe)	max.0.001%

JIS K 8701-1994
*:Additional test performed by Wako but not required by JIS
This product is shot