

Code No.12-0046

SPECIFICATION

Lead(II) Acetate Trihydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White crystals
Assay [(CH ₃ COO) ₂ Pb·3H ₂ O]	min.99.5%(mass/mass)
Solubility in dil.acetic acid	to pass test
Chloride (Cl)	max.3ppm(mass/mass)
Nitrate	to pass test
Sodium (Na)	max.0.005%(mass/mass)
Potassium (K)	max.0.005%(mass/mass)
Copper (Cu)	max.5ppm(mass/mass)
Magnesium (Mg)	max.0.005%(mass/mass)
Calcium (Ca)	max.0.005%(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)

JIS K 8374:2007

*:Additional test performed by Wako but not required by JIS