

## SPECIFICATION

## Lead(II) Nitrate

## Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay [Pb(NO <sub>3</sub> ) <sub>2</sub> ]	min.99.5%(mass/mass)
Solubility in water	to pass test
Chloride (Cl)	max.0.001%(mass/mass)
Sodium (Na)	max.0.005%(mass/mass)
Potassium (K)	max.0.005%(mass/mass)
Copper (Cu)	max.0.001%(mass/mass)
Calcium (Ca)	max.0.005%(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)

JIS K 8563 : 2018

(\*:Additional test performed by Wako but not required by JIS)