

Code No.12-0070

## SPECIFICATION

### Dilead(II) Lead(IV) Oxide

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	Reddish yellow ~ red, powder
Solubility in nitric acid-hydrogen peroxide	to pass test
Loss on drying at 105°C	max.0.05%
Chloride (Cl)	max.0.005%
Nitrate	to pass test(NO3:max.about 0.005%)
Manganese (Mn)	max.0.005%
Assay	min.97.0%