

Code No.12-0085

SPECIFICATION

L-Leucine

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	White, crystals ~ crystalline powder
Assay (C ₆ H ₁₃ NO ₂)	min.99.0%(mass/mass)
Solubility in dil.hydrochloric acid	to pass test
pH of a 10g/L solution at 25°C	5.5~7.0
Specific rotation [α] _{D20}	+14.7~+16.0°
Loss on drying at 105°C	max.0.2%(mass/mass)
Residue after ignition (as sulfate)	max.0.05%(mass/mass)
Chloride (Cl)	max.0.01%(mass/mass)
Sulfate (SO ₄)	max.0.02%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Arsenic (As)	max.1ppm(mass/mass)
Iron (Fe)	max.0.001%(mass/mass)
Ammonium (NH ₄)	max.0.02%(mass/mass)
Foreign amino acids	to pass test

JIS K 9054 : 2020(*:Additional test performed by Wako but not required by JIS)