

Code No.12-0153

SPECIFICATION

Lithium Acetate Dihydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder or mass to pass test
Solubility in water	
pH of a 50g/l solution at 25°C	7.0~9.4
Chloride (Cl)	max.0.005%
Heavy metals (as Pb)	max.0.001%
Sodium (Na)	max.0.02%
Potassium (K)	max.0.005%
Iron (Fe)	max.0.001%
Assay	98.0~102.0%