

# SPECIFICATION

## Lithium Tetraborate, Anhydrous

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, powder or mass
Bulk density	min.0.5g/mL
Loss on ignition	max.1.0%
Chloride (Cl)	max.0.01%
Sulfate (SO <sub>4</sub> )	max.0.02%
Magnesium (Mg)	max.0.005%
Calcium (Ca)	max.0.005%
Manganese (Mn)	max.0.005%
Iron (Fe)	max.0.005%
Assay	min.98.0%