

SPECIFICATION

Lanthanum Chloride Heptahydrate

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White crystals
Solubility in water	to pass test
Nitrate (NO ₃)	max.0.01%
Sulfate (SO ₄)	max.0.01%
Sodium (Na)	max.0.01%
Potassium (K)	max.0.005%
Copper (Cu)	max.0.001%
Magnesium (Mg)	max.0.001%
Calcium (Ca)	max.0.01%
Cerium (Ce)	max.0.05%
Praseodymium (Pr)	max.0.05%
Neodymium (Nd)	max.0.05%
Lead (Pb)	max.0.002%
Iron (Fe)	max.0.001%
Assay	min.97.0%