

Code No.12-0687

SPECIFICATION

Lysyl EndopeptidaseR, recombinant, Biopharmaceutical Analysis Grade
for Cell Culture

| REQUIREMENT | SPECIFICATION |
|---------------------------------|-------------------|
| Appearance | Lyophilisate |
| Bacterial endotoxins | less than 20EU/mg |
| Residual DNA | max.10ng/mg |
| Host cell-derived protein | max.1 μ g/mg |
| Electrophoresis test (SDS-PAGE) | to pass test |
| Activity | min.2.0AU/mg |

Source:E. coli expressed Lysyl EndopeptidaseR