

SPECIFICATION

Manganese(II) Nitrate Hexahydrate

Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	to pass test
Assay[Mn(NO ₃) ₂ ·6H ₂ O]	min.98.0%(mass/mass)
Solubility in water	to pass test
pH of a 50g/l solution at 25°C	3.0~4.5
Chloride (Cl)	max.0.001%(mass/mass)
Sulfate (SO ₄)	max.0.002%(mass/mass)
Sodium (Na)	max.0.005%(mass/mass)
Potassium (K)	max.0.002%(mass/mass)
Copper (Cu)	max.5ppm(mass/mass)
Magnesium (Mg)	max.0.005%(mass/mass)
Calcium (Ca)	max.0.005%(mass/mass)
Zinc (Zn)	max.0.001%(mass/mass)
Lead (Pb)	max.0.001%(mass/mass)
Iron (Fe)	max.5ppm(mass/mass)
Nickel (Ni)	max.0.002%(mass/mass)

JIS K 8568:2011(*:Additional test performed by Wako but not required by JIS)