

Code No.13-0102

SPECIFICATION

Mercury
Guaranteed Reagent

REQUIREMENT	SPECIFICATION
* Appearance	to pass test
Assay	min.99.5%
Solubility in nitric acid	to pass test
Residue after evaporation	max.5ppm
Copper (Cu)	max.0.5ppm
Zinc (Zn)	max.0.5ppm
Lead (Pb)	max.0.5ppm
Iron (Fe)	max.0.5ppm

JIS K 8572-1992

*:Additional test performed by Wako but not required by JIS