

Code No.13-0160

## SPECIFICATION

## L-Methionine

Wako Special Grade

REQUIREMENT	SPECIFICATION
Appearance	White, crystals ~ crystalline powder
Solubility in water	to pass test
Specific rotation [ $\alpha$ ]D20(c=8,HCl) (after drying)	+23.7~+24.7°
Loss on drying at 105°C	max.0.2%
Chloride (Cl)	max.0.01%
Sulfate (SO <sub>4</sub> )	max.0.02%
Heavy metals (as Pb)	max.0.002%
Arsenic (As)	max.1 ppm
Iron (Fe)	max.0.002%
Ammonium (NH <sub>4</sub> )	max.0.02%
Thin layer chromatography test	to pass test
Assay (after drying)	min.99.0%